

## Cameras Currently Supported by Aurox Visionary Software (Version 4.2.9)

Camera model	Features supported in Visionary 4.2.9
<b>Hamamatsu Orca Flash 4 V3</b>	<p>USB 3.0 or Camera Link Connection available.</p> <p>Up to 50 fps acquisition with USB connection or up to 100 fps acquisition with camera link connection (for use with Clarity HS).</p> <p>Binning and frame averaging supported.</p>
<b>PCO Edge 5.5 sCMOS</b>	<p>USB 3.0 connection.</p> <p>Up to 37 fps acquisition with USB connection.</p> <p>Binning and frame averaging supported.</p>
<b>Andor Zyla 5.5 sCMOS</b>	<p>USB 3.0 connection.</p> <p>Up to 37 fps acquisition with USB connection.</p> <p>Binning and frame averaging supported.</p>
<b>Photometrics 95PrimeB</b>	<p>USB 3.0 Connection.</p> <p>Up to 40 fps acquisition with USB connection.</p> <p>Binning and frame averaging supported.</p>
<b>Hamamatsu Orca Flash 4 V2 (USB 3.0)</b>	Some integration and testing complete.
<b>Hamamatsu Orca Flash 4 LT</b>	Some integration and testing complete.
<b>PCO Panda</b>	Some integration and testing complete.

If other cameras are required, please contact Aurox to investigate the possibility of integration