

## Light Sources Currently Supported by Aurox Visionary Software (Version 4.2.9)

Light source model	Features supported in Visionary 4.2.9
<b>CoolLED pE300<sup>ultra</sup></b> (Preferred light source of choice)	On/off control as a white light source. Global intensity control Coupled directly to Clarity LFC using Zeiss adapter.
<b>CoolLED pE4000</b>	On/off control as a white light source. Global intensity control Coupled directly to Clarity LFC using Zeiss adapter.
<b>Excelitas X-cite 120 PC-Q</b>	On/off control during acquisition as a white light source and global intensity control.
<b>Andor AMH200</b>	On and off control during acquisition as a white light source.
<b>Lumencor Sola</b>	On/off control only .
<b>Lumencor Spectra-X</b>	On/off control and intensity control.
<b>Prior Lumen 200</b>	On/off control as white light source.
<b>Excelitas X-cite 120 LED</b>	On/off control during acquisition as a white light source and global intensity control.

If other light sources are required, please contact Aurox to investigate the possibility of integration